

Oral Fluid Drug Screening Now Available!

Fast and Simple Collection Procedure

FDA approved Oral Fluid Collection Device purchased separately.

See package insert for complete instructions.

WARNING: Do not drink the blue fluid in the vial.

Step 1: Swab – Swab for 5 minutes.



Step 2: Snap – Insert pad in vial. Snap wand at scoring.



Step 3: Seal – Recap and apply seal. Send to STERLING for Drug Screening!



Customize from 4 to 11 Drugs:

Amphetamine
Methamphetamine
Benzodiazepine
Barbiturate
Cocaine
Cannabinoids (marijuana)
Methadone
Opiates
Oxycodone
Propoxyphene
Phencyclidine (PCP)

Confirmation available By Request Only,
or Automatically, at a separate fee.

***Negative AND Positive
Screen results within 24
hours from receipt at lab!***

Contact our Sales Department for more information!

1-800-442-0438 or sales@regtox.com